

Injection Date : 9/23/97 11:04:00 AM Seq. Line : 1
Sample Name : STD1 Vial : 1
Acq. Operator : martinj Inj : 1
Inj Volume : 1 µl
Acq. Method : C:\HPCHEM\1\METHODS\H2O.M
Last changed : 9/23/97 11:00:15 AM by martinj
Analysis Method : C:\HPCHEM\1\METHODS\NIC.M
Last changed : 10/3/97 10:41:00 AM
(modified after loading)

Water Analysis - Routine Nicotine and Water Analysis

-----
Internal Standard Report (Sample Amount is 0!)

Sorted By : Signal
Calib. Data Modified : 10/3/97 10:35:20 AM
Multiplier : 1.0000
Dilution : 1.0000
Uncalibrated Peaks : not reported

Sample ISTD Information:

ISTD ISTD Name
[mg/cigt]

1 1.00000 Carvone

Signal 1: FID2 B,
Results obtained with enhanced integrator!

RetTime [min]	Type	Area counts*s	Amt/Area ratio	Amount [mg/cigt]	Grp	Name
2.819	VV X+I	2.93162e4	1.00000	1.00000		Carvone
3.683	VV T+	2923.34448	2.91928e-1	2.91104e-2		Nic

Totals without ISTD(s) : 2.91104e-2

1 Warnings or Errors :

Warning : Sample amount is zero. Absolute amounts calculated